


marClamp[®] CUT IQ



marClamp[®] CUT IQ

Unbeatable for Efficiency and Economy



marClamp® CUT IQ

The new combination instrument for open surgical vessel sealing

maxium®

+ *marClamp*® CUT IQ

= SealSafe® IQ

... the more intelligent sealing system!



All the benefits at a glance

The *marClamp*® CUT IQ instruments, in combination with the maxium® electrosurgical unit, are used for bipolar vessel sealing, followed by tissue separation, **with no need to change instruments!**

Convincingly simple due to “Plug & Play” function

The maxium® HF generator recognizes the connected instrument thanks to innovative IQ technology and sets the required parameters immediately.

As soon as the tissue is gripped, the switch integrated into the instrument automatically triggers the special SealSafe® IQ sealing current so no additional foot switch is required.

When sealing has been performed, the current is stopped and an acoustic signal sounds.

Economically unbeatable

The ingenious part about *marClamp*® CUT IQ is the very low cost per application, which is the result of complete reusability and the minimal proportion of disposables accounted for by the sterile blade.

- Reliable dissection, sealing and separation with just one instrument
- Reproducible, swift sealing process
- Shorter operating times
- Minimal lateral heat spread
- Substantial cost saving due to reusability and low proportion accounted for by disposables
- Optimal cutting results due to sterile disposable blade
- Manual activation
- “Plug & Play” function
- Easy handling and reprocessing
- Autoclavable at 134°C
- Sealing length: 17 mm
- Cutting length: 15 mm




“Reuse makes economic sense. We value safety, reliability, and cost-effectiveness!”

Prof. Gerald Gitsch
University Women's Hospital Freiburg
Medical Director



Cost comparison



Applications

Gynaecology

- Lymphadenectomy
- Adhesiolysis
- Cystectomy

Visceral surgery

- Colectomy
- Gastrectomy
- Splenectomy
- Adrenalectomy
- Nephrectomy

ENT

- Thyroidectomy
- Neck dissection
- Laryngectomy

Urology

- Cystectomy
- Hydrocele
- Partial cystectomy

Article overview

Art. no.	Description
80-632-00-04	Blade for marClamp® CUT IQ, sterile, pack of 10
80-632-18-04	marClamp® CUT IQ set, 18 cm
80-630-58-04	Cleaning tool
80-632-01-04	Set of cleaning brushes, 3 pcs.

KLS Martin Group

Karl Leibinger Medizintechnik GmbH & Co. KG
78570 Mühlheim · Germany
Tel. +49 7463 838-0
info@klsmartin.com

KLS Martin GmbH + Co. KG
79224 Umkirch · Germany
Tel. +49 7665 9802-0
info@klsmartin.com

Stuckenbrock Medizintechnik GmbH
78532 Tuttlingen · Germany
Tel. +49 7461 165880
verwaltung@stuckenbrock.de

Rudolf Buck GmbH
78570 Mühlheim · Germany
Tel. +49 7463 99516-30
info@klsmartin.com

KLS Martin France SARL
68200 Mulhouse · France
Tel. +33 3 8951 3150
france@klsmartin.com

Martin Italia S.r.l.
20864 Agrate Brianza (MB) · Italy
Tel. +39 039 605 6731
italia@klsmartin.com

Martin Nederland/Marned B.V.
1271 AG Huizen · The Netherlands
Tel. +31 35 523 4538
nederland@klsmartin.com

KLS Martin UK Ltd.
Reading RG1 3EU · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

Nippon Martin K.K.
Osaka 541-0046 · Japan
Tel. +81 6 62 28 9075
nippon@klsmartin.com

KLS Martin L.P.
Jacksonville, FL 32246 · USA
Tel. +1 904 641 7746
usa@klsmartin.com

KLS Martin do Brasil Ltda.
CEP 04.531-011 São Paulo · Brazil
Tel.: +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Australia Pty Limited
Artarmon NSW 2064 · Australia
Tel.: +61 2 9439 5316
australia@klsmartin.com

KLS Martin Malaysia Sdn. Bhd.
10200 Penang · Malaysia
Tel. +604 263 2566
malaysia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
121471 Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
201203 Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Postfach 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

